

Product Information

CTA: 146189

roboTest L (Linear) robotic testing system for plastic



roboTest L robotic testing system for plastics

Range of application

The system is used for the following fully automatic tests:

• tensile and flexure tests on plastics specimens (e.g. to ISO527-2 (1A, 1B, 1BA), ISO 178 or ASTM D790)

System configuration

- 5kN or 250kN materials testing machine with symmetrically closing pneumatic or hydraulic grips and optional extensometer
- roboTest L automatic specimen feed system
- Cross-section measuring device (optional)
- Industrial controller with testXpert testing software and autoEdition3 automation software

Benefits

- ZwickRoell has over 35 years of experience and expertise, gained while supplying more than 700 automated testing systems worldwide.
- Operator influences (hand temperature/moisture, offcenter or angled specimen insertion, etc.) are eliminated for high test-result reproducibility.
- Qualified laboratory staff are relieved of routine activities, making them available for more complex tasks.
- The machine can be used during idle times (lunch breaks and night shifts), which increases capacity and produces faster results.

- The testing system reduces the testing costs per specimen and typically pays for itself within one to two years.
- Manual tests can be performed whenever required the specimen feeder simply slides out of the way for access.

Test sequence

- The operator fills the removable specimen magazine during specimen preparation or directly on the system. The magazine is then placed on the magazine table.
- Specimen data (identification number, width, thickness etc.) are entered on the PC. This step can be omitted when using barcodes.
- Once the system has been started, specimen feed, tensile testing and disposal of specimen remains are performed automatically.
- Once all the specimens in the magazine have been tested it can be refilled or replaced with a magazine insert which has already been filled. Alternatively magazines can be 'topped up' with specimens during the automatic sequence.



Product Information

roboTest L (Linear) robotic testing system for plastic

Technical data

Туре	roboTest L		
Mechanics			
Attachment	Can be docked to the load frame	Can be docked to the load frame	
Dimensions (H x W x D)			
roboTest L, metals	1600 ¹⁾ x 1330 ²⁾ x 1500 ³⁾	mm	
roboTest L, plastics	1200 ¹⁾ x 1230 ²⁾ x 1600 ³⁾	mm	
Weight (without specimens, with cross-section meas	suring device and magazine)		
roboTest L, metals	250	kg	
roboTest L, plastics	200	kg	
Power specifications			
Nominal supply voltage	230	V AC	
Power consumption	200	Procedure Instruc- tion	
Power frequency	50/60	Hz	
Compressed air, filtered, not oiled	5 6	bar	
Compressed air consumption, approx.	30	l/h	
Controller			
Automation	autoEdition3		
Peripheral connection	Profinet		

1) Dependent on the testing machine

2) Incl. movable magazine table; additional travel distance: 990 mm

3) With linear axis

Value	
Plastics	
tensile tests, flexure tests	
dumbbell, strip	
pincer or suction gripper	
max. 600	
rigid, non-tacky	
max. 0.1	kg
max. 260	mm
6 to 25	mm
-	mm
max. 15	mm
	Plastics tensile tests, flexure tests dumbbell, strip pincer or suction gripper max. 600 rigid, non-tacky max. 0.1 max. 260 6 to 25

Description	ArticleNumber
roboTest L automatic specimen feed system with pincer gripper	3008305
For tensile tests	
 Pneumatic specimen gripper with rotation unit Magazine table with centering pins for 4 magazine inserts 	

All data at ambient temperature.

Zwick Roell

Product Information

roboTest L (Linear) robotic testing system for plastic

Description	ArticleNumber
Technical data • Gripper opening width: max. 20 mm • Specimen thickness 0,1 to 18 mm • Air consumption: approx. 30 l/h	
 roboTest L automatic specimen feed system with suction gripper For tensile and flexure tests Magazine table with centering pins for 48 specimen-stacks in 4 rows Technical data Specimen thickness 1 to 6 mm Air consumption approx. 30 l/h 	1090735
Options	

- Specimen identification
- Cross-section measurement No. measurements/specimen: 1/3
- Specimen disposal
- Good/bad sorter
- Data exchange: Higher-level computer system (e.g. LIMS) via upload/download of ASCII filesor ODBC
- Visual status display: 3-aspect lights (running, refill magazine/finished, error)