

Product Information

Roller Grips, Type 8564 Fmax 100, Type 8566 Fmax 250 kN





Type 8566 roller grips, Fmax 250 kN

Applications

- Specimen material: textiles
- Specimen shape: strips (belt, strap)
- Type of loading: tensile

Function description

These roller grips are used for testing tensile strip specimens.

In tests with these grips optical extensometers are used for accurate strain measurement.

The self-clamping action of the grips is generated by multiple looping of the specimen and the application of tensile force.

General view: Type 8566 roller grips, Fmax 100 kN and Fmax 250 kN

The specimen is inserted into a split roller from the front. The looping occurs with the rotation of the roller.

The specimen grip locking bars are monitored for safety in conjunction with the machine electronics.

Advantages and features

- Ergonomic design for fast, easy specimen insertion
- Low overall height of the specimen grips allows specimens with large elongations to be tested.
- Suitable for tests in temperature chambers
- The grips are suitable for clamping-sensitive specimens.



Product Information

Roller Grips, Type 8564 Fmax 100, Type 8566 Fmax 250 kN

Technical data

313748 ¹⁾	313840 ¹⁾	
8564	8566	
100	250	kN
210	270	mm
260	400	mm
170	230	mm
233	305	mm
100	130	mm
110	220	mm
Steel, pr ²⁾ 1 mm	Steel, pr ²⁾ 1.5 mm	
8	10	mm
110	220	mm
1000	1000	mm
24	63	kg
Ø 60	Ø 60	mm
+10 +35	+10 +35	°C
2	2	pieces
	313748¹⁾ 8564 100 210 260 170 233 100 110 Steel, pr ²⁾ 1 mm 8 110 1000 24 24 24 0 60 +10 +35 2	3137481) 3138401) 8564 8566 100 250 210 270 260 400 170 230 233 305 110 220 Steel, pr ²) 1 mm Steel, pr ²) 1.5 mm 8 10 110 220 Steel, pr ²) 1 mm Steel, pr ²) 1.5 mm 100 1000 110 20 110 210 110 100 110 20 110 20 110 20 110 20 110 20 110 20 110 20 110 20 1000 1000 1000 0 100 0 100 1000 24 63 100 10 100 10 100 10 10

1) Only in connection with Xforce K load cell

2) Abbreviations: Pr = pyramid grid