

## **Product Information**

TorsionLine TL 020 / TL 200 / TL 500 Torsion Testing Machines



TorsionLine torsion testing machine with testControl II

### **Range of application**

TorsionLine TL020, TL200 and TL500 torsion testing machines are used to test torsional influence on materials and components. The torsion testing axis is horizontally orientated.

### Main areas of use


- Main areas of use
- Plastics

CTA: 118031 11757

- Metallic materials, e.g. reverse torsion test on wire, EN 2002-13
- Medical engineering products, e.g. metal bonescrews, ISO 6475
- Fasteners and connectors
- Bolts/screws in general
- Composites
- Drive shafts
- Motor/engine bearings

### **Advantages and Features**

- Maximum stiffness ensures precise rotation angle measurement over the entire torque range.
- High resolution of rotation angle measurement



Drawing: TorsionLine with testControl II

- Use of maintenance-free AC servo drives
- Operation with standard PC (no additional interface card required) and testXpert III testing software
- Designed for production, quality assurance, and research and development
- Maximum flexibility in handling and operation through state-of-the-art software control
- Safety device with electromechanical guard locking (CE compliant)
- ZwickRoell comprehensive range of optional accessories
- Mechanical overload protection for small torque transducers during specimen installation
- Large test speed range for high specimen throughput
- Optional static weight loading unit (tensile and compression) for application of axial forces up to max. 500 N

# Zwick Roell

## **Product Information**

TorsionLine TL 020 / TL 200 / TL 500 Torsion Testing Machines

Туре	TL 020	TL 200	TL 500	
Item No.	1026875	1026878	1026879	
Load frame				
Torque (left/right)	20	200	500	Nm
Total height with safety device	1775	1775	1775	mm
Total width with safety cover	1475	1475	1475	mm
Total depth	650	650	650	mm
Height of test axis	Adjustable	Adjustable	Adjustable	
Safety device	Electrically inter- locked	Electrically inter- locked	Electrically inter- locked	
Anti-rotation device of torque transducer	Manually operated, electrically monitored	Manually operated, electrically monitored	Manually operated, electrically monitored	
Weight with electronics, safety cover, tor- que transducer	185	188	188	kg
Grip to grip separation with 4-jaw chuck, max.	491	491	491	mm
Grip-to-grip separation with universal flanges (without specimen grip), max.	680	680	680	mm
Diameter of specimens/specimen hold- ers, max.	200	200	200	mm
Finish	RAL 7011 and RAL 7038	RAL 7011 and RAL 7038	RAL 7011 and RAL 7038	
Ambient temperature	+10 +35	+10 +35	+10 +35	°C
Relative humidity (non-condensing)	20 90	20 90	20 90	%
Noise level	70	70	70	dB(A)
Torsion inherent stiffness	> 1000	> 1000	> 1000	Nm/°
Axial frictional force when driving the car- riage (without torsion load), approx.	4	4	4	Ν
Drive				
Rotational speed				
In turns	0.0005 180	0.0005 90	0.0005 50	rpm
In degrees	0.072 64800	0.036 32400	0.036 18000	°/min
Angle resolution of the drive	3.29	1.89	0.92	Rotational second/ minute
Positioning repeatability (without revers- ing directions)				
At 1000 °/min	≤20	≤10	≤10	Rotational seconds
Play at zero torque crossing	< 1	< 1	< 1	Rotational minutes
Power specifications				
Electrical connection	230V, 1 Ph/N/PE	400V, 3 Ph/N/PE	400V, 3 Ph/N/PE	
Power consumption	2.2	7	7	kVA
Frequency	50/60	50/60	50/60	Hz



## **Product Information**

## TorsionLine TL 020 / TL 200 / TL 500 Torsion Testing Machines

#### **Torque transducer**

Anti-rotation device included in delivery

For TorsionLine	Nominal torque M <sub>nom</sub> [Nm]	Measurement range/ measurement travel	Item No.
		Accuracy Class 1 <sup>1)</sup> [Nm]	
TL 020, TL 200	2	≥0.02	623682
TL 020, TL 200	5	≥0.05	050251
TL 020, TL 200	20	≥0.2	623681
TL 200	200	≥2	623680
TL 500	500	≥5	631051

1) Based on ISO 7500-1