


Paper strip cutter

CTA: 46349



Applications

All paper strip cutters are supplied with hold-down clamp and a set of spare cutting-wheels. They are suitable for papers with a grammage from 30 g/m² to 350 g/m².

Paper strip cutter

Technical data

Туре	Paper strip cutter	Paper strip cutter	Paper strip cutter	Paper strip cutter	
Item No.	377820	320412	363795	320422	
For strip width	12.7	15	25.4	50	mm
Length of cut, max.	350	350	350	350	mm
Ambient temperature	+10 +35	+10 +35	+10 +35	+10 +35	°C

Optional accessories

Description	ArticleNumber
Replacement counter-blade, to fit paper strip cutter with strip width 12.7 mm	377821
Replacement counter-blade, to fit paper strip cutter with strip width 15 mm	320414
Replacement counter-blade, to fit paper strip cutter with strip width 50 mm	320424